

## Setting Up a Tissue Culture Lab:

---

Any laboratory, in which tissue culture techniques are performed, regardless of the specific purpose, must contain a number of basic facilities. These usually include the following:

- A general washing area
- A media preparation, sterilization, and storage area
- An aseptic transfer area
- Environmentally controlled incubators or culture rooms
- An observation/data collection area.

### Washing Area

The washing area should contain large sinks, some lead-lined to resist acids and alkalis, draining boards, and racks, and have access to demineralized water, distilled water, and double-distilled water. Space for drying ovens or racks, automated dishwashers, acid baths, pipette washers and driers, and storage cabinets should also be available in the washing area.

### Media Preparation Area

The media preparation area should have ample storage space for the chemicals, culture vessels and closures, and glassware required for media preparation and dispensing. Bench space for hot plates/stirrers, pH meters, balances, water baths, and media-dispensing equipment should be available. Other necessary equipment may include air and vacuum sources, distilled and double-distilled water, Bunsen burners with a gas source, refrigerators and freezers for storing stock solutions and chemicals, a microwave or a convection oven, and an autoclave or domestic pressure cooker for sterilizing media, glassware, and instruments.

In preparing culture media, analytical grade chemicals should be used and good weighing habits practiced. To insure accuracy, and exact step-by-step routine should be developed for media preparation and a complete checklist required of all media preparers even for the simplest media.

The water used in preparing media must be of the utmost purity and highest quality. Tap water is unsuitable because it may contain cations (ammonium, calcium, iron, magnesium, sodium, etc.), anions (bicarbonates, chlorides, fluorides, phosphates, etc.), microorganisms (algae, fungi, bacterial), gases (oxygen, carbon dioxide, nitrogen), and particulate matter (silt, oils, organic matter, etc.) Water used for plant tissue culture should meet, at a minimum, the standards for type II reagent grade water, i.e., be free of pyrogens, gases, and organic matter and have an electrical conductivity less than 1.0  $\mu\text{mho/cm}$ .

The most common and preferred method of purifying water to type II standards is a deionization treatment followed by one or two glass distillations. The deionization treatment removes most ionic impurities, and the distillation process removes large organic molecules, microorganisms, and pyrogens. Three other methods that will produce type II purity water are absorption filtration, which uses activated carbon to remove organic contaminants and free chlorine; membrane filtration, which removes particulate matter and most bacterial contamination; and reverse osmosis, which removes approximately 9% of the bacterial, organic, and particulate matter as well as about 90% of the ionized impurities.

### Transfer Area

Under very clean and dry conditions, tissue culture techniques can be successfully performed on an open laboratory bench. However, it is advisable that a laminar flow hood or sterile transfer room be utilized for making transfers. Within the transfer area there should be a source of electricity, gas, compressed air, and vacuum.

The most desirable arrangement is a small dust-free room equipped with an overhead ultraviolet light and a positive-pressure ventilation unit. The ventilation should be equipped with a high-efficiency particulate air (HEPA) filter. A 0.3- $\mu\text{m}$  HEPA filter of 99.97-99.99% efficiency works well. All surfaces in the room should be designed and constructed in such a manner that dust and microorganisms do not accumulate and the surfaces can be thoroughly cleaned and disinfected. A room of such design is particularly useful if large numbers of cultures are being manipulated or large pieces of equipment are being utilized.

Another type of transfer area is a laminar flow hood. Air is forced into the unit through a dust filter then passed through a HEPA filter. The air is then either directed downward (vertical flow unit) or outward (horizontal flow unit) over the working surface. The constant flow of bacteria-free filtered air prevents nonfiltered air and particulate matter from settling on the working surface.

The simplest type of transfer area suitable for tissue culture work is an enclosed plastic box commonly called a glove box. This type of culture hood is sterilized by an ultraviolet light and wiped down periodically with 95% ethyl alcohol when in use. This type of unit is used when relatively few transfers are required.

### Culture Room

All types of tissue cultures should be incubated under conditions of well-controlled temperature, humidity, air circulation, and light quality and duration. These environmental factors may influence the growth and differentiation process directly during culture or indirectly by affecting their response in subsequent generations. Protoplast cultures, low-density cell suspension cultures, and anther cultures are particularly sensitive to environmental cultural condition.

Typically, the culture room for growth of plant tissue cultures should have a temperature between 15° and 30° C, with a temperature fluctuation of less than  $\pm 0.5^\circ\text{C}$ ; however, a wider range in temperature may be required for specific experiments. It is also recommended that the room have an alarm system to indicate when the temperature has reached preset high or low temperature limits, as well as continuous temperature recorder to monitor temperature fluctuations. The temperature should be constant throughout the entire culture room (i.e., no hot or cold spots). The culture room should have enough fluorescent lighting to reach the 10,000 lux; the lighting should be adjustable in terms of quantity and photoperiod duration. Both light and temperature should be programmable for a 24-hr period. The culture room should have fairly uniform forced-air ventilation, and a humidity range of 20-98% controllable to  $\pm 3$  percent. Many incubators, large growth chambers, and walk-in environmental chambers meet these specifications.

## BASIC LABORATORY EQUIPMENT

Many tissue culture techniques require similar basic laboratory equipment. The following items are commonly found in a laboratory for *in vitro* propagation of plant materials:

### Items Required to Start a Plant Cell Culture Laboratory

QTY	ITEM DESCRIPTION	PHYTOTECH CATALOGUE NO.	APPROX. COST	ITEM FUNCTION
1	Water Purification System; water should have a resistivity of at least 200,000 ohms-cm and a conductivity 5.0 micromhos/cm	NA		Purification of water for media preparation
1	Electronic Balance (0.01 g readability; 200 g minimum capacity)	B 933	295.00	Measuring out biochemicals and media
1	pH meter (range 0-14 +/- 0.01; automatic temperature compensation 0-60 <sup>0</sup> C; one or two point calibration)	P 976	74.95	Measurement and adjustment of media pH
1	Hot Plate/Stirrer (7" x 7" ceramic top; variable heating range from ambient to 400 <sup>0</sup> C; variable stirring speed from 50-150 rpm; chemically resistant)	H 926	295.00	Mixing & heating media and stock
1	Refrigerator/freezer; capable of maintaining a refrigerator temperature of 0-5 <sup>0</sup> C with a freezer temperature of approximately -20 <sup>0</sup> C	NA	750.00	Storage of stock solutions, media, hormones
1	Laminar Flow Transfer Hood; incoming air should be HEPA filtered to remove 99.99% of particles larger than 0.3µm; should meet or exceed the Class 100 Clean Standard 209D; maintain a flow of 90 fpm +/- 20% at static pressures of 0.6-1.2"	NA	3500.00	Provide a sterile atmosphere to transfer cultures
	4 liter Isopropyl alcohol	NA	10.00	Used to sterilize instruments and work areas
1 roll	Aluminum foil, heavy duty; (18" x100 ft roll)	NA	2.50	Used to wrap instruments prior to sterilization, cover vessels
12	Beakers, 250 mL	B 910	7.75	Mixing solutions
12	Beakers, 1000 mL	B 931	13.75	Mixing solutions
6	Beakers, 2000 mL	B 939	32.25	Mixing solutions and media
6	Beakers, 4000 mL	B 960	46.25	Mixing solutions and media
4 ea	Bottle, Water; pre-labeled for use in dispensing water; 500 mL capacity	B 974	6.00	Rinsing instruments, beakers, transplants from tissue culture
4 ea	Bottle, Isopropyl Alcohol; pre-labeled for use in dispensing isopropyl alcohol; 500 mL capacity	B 987	6.00	Rinsing sterile hood work surfaces
1 case	Bottle, 100 mL; Type 1 borosilicate glass with volume graduations; supplied with 33-340 black polypropylene cap with rubber liner; non-sterile		24.00	Storage of stock solutions, sterile distilled water, media
1 case	Bottle, 500 mL; Type 1 borosilicate glass with volume graduations; supplied with 33-340 black polypropylene cap with rubber liner; non-sterile		50.00	Storage of stock solutions, sterile distilled water, media
2	Brush, flask or bottle			Cleaning glassware

## Items Required to Start a Plant Cell Culture Laboratory

QTY	ITEM DESCRIPTION	PHYTOTECH CATALOGUE NO.	APPROX. COST	ITEM FUNCTION
1 case	Culture tubes, 25 x 150 mm, borosilicate glass; 500 tubes/case	C 930	146.50	Starting cultures in Stage I
1 case	Culture tube racks; holds 40, 25 mm culture tubes; withstands temperatures up to 121 <sup>0</sup> C	C 908	90.50	Holding culture tubes
500	Closure, for 25 mm culture tubes, 500 each	C 925	57.50	Sealing culture tubes
1 case	Culture vessel, baby food jar; glass culture vessel; autoclaveable; uses Magenta B Cap (C 903) as closure; 110 mm height; 100/case	C 904	42.50	Culture vessel for Maintaining plant cultures
1 case	Culture vessel, baby food jar; glass culture vessel; autoclavable; uses Magenta B Cap (C 903) as closure; 72 mm height; 100/case	C 900	32.50	Culture vessel for Maintaining plant cultures
1 case	Magenta B Cap; autoclavable closure for baby food jars; fits both C 904 and C 900; clear polypropylene closure; 100/case	C 903	215.50	Closure for baby food culture Vessel
2 cases	Culture vessel; autoclavable culture vessel and lid made from clear polypropylene; round vessel measures; 250/case	C 913		Culture vessel for maintaining plant cultures
1 gal	Detergent	NA	5.75	Cleaning glassware
1 case	Culture dishes, disposable, sterile, 100 x 25 mm	D 940	173.50	For Stage I cultures sterile surface for cutting explants
1 gal	Disinfectant, commercial brand	NA	2.00	Disinfects explants
1	Erlenmeyer flask, wide mouth, 1000 mL	F 985	18.50	Mixing media
1	Erlenmeyer flask, wide mouth, 2000 mL	F 986	25.50	Mixing media
1	Erlenmeyer flask, 4000 mL			Mixing media
1	Erlenmeyer flask, 6000 mL			Mixing media
1 case	Filtration system, vacuum; disposable, plastic, sterile system designed for filtration of fluids needed in tissue culture; polystyrene, screw-cap base for sterile storage; 200 mL; 47 mm diameter/0.22 µm pore size nylon membrane; 12 units/case	NA	63.00	Sterilization of heat liable stock solutions
3	Forceps, dressing; 10" length, serrated, stainless steel	F 952	30.75	Transferring tissue
3	Forceps, bayonet; stainless steel	F 957	71.85	Transferring tissue
5	Forceps; stainless steel, very fine point, 6" length			Transferring tissue
1	Graduated cylinders, glass or polypropylene, 10 mL	C 962	3.50	Preparing stock solutions
1	Graduated cylinders, glass or polypropylene, 100 mL	C 964	6.25	Preparing stock solutions
1	Graduated cylinders, glass or polypropylene, 1000 mL	C 968	20.25	Preparing stock solutions
1 pkg	Lab markers, assorted colors (10)	NA		Labeling cultures

## Items Required to Start a Plant Cell Culture Laboratory

QTY	ITEM DESCRIPTION	PHYTOTECH CATALOGUE NO.	APPROX. COST	ITEM FUNCTION
1 roll	Parafilm (4" x 250 ft)	S 911	36.50	Wrapping culture Closure
100	Pipets, 1 mL, graduated, sterile, disposable, individually wrapped	P990	14.50	Measuring out stock solutions
100	Pipets, 5 mL, graduated, sterile, disposable, individually wrapped	P993	26.25	Measuring out stock solutions
100 ea	Pipets, 10 mL, graduated, sterile, disposable, individually wrapped	P994	24.95	Measuring out stock solutions
100 ea	Pipets, 25 mL, graduated, sterile, disposable, individually wrapped	P995	58.25	Measuring out stock solutions
1	Pipet pump, electric	NA	190.00	For safely measuring liquids
1 pkg	Gloves, hot or autoclave; provides protection up to 350 <sup>0</sup> F	NA	6.00	Safely removing hot items from autoclave
2 ea	Scalpel handle; No 3: 5" length; stainless steel	S963	69.75	Cutting explants
2 ea	Scalpel handle; No 3L: 8" length; stainless steel	S973	83.75	Cutting explants
1 box	Scalpel blades; No 10; stainless steel; individually wrapped; sterile; 100/box	S970	56.25	Cutting explants
1	Scoop, large plastic 5-3/4" x 9" bowl	NA	11.00	Measuring large volumes of biochemical's
1 pkg	Scoop, medium, 30 cm length, 2/pkg	NA	34.65	Measuring small to medium amounts of biochemical's
2 ea	Spatula, stainless steel blade, wooden handle, 4" wide x 19/32" thick	NA	11.00	Measuring small to medium amounts of biochemical's
1 pkg	Spatula, micro; 8" length, flat nickel stainless steel	NA	25.25	Measuring small to medium amounts rounded ends, 2" end length; of biochemical's
1 ea	Sterilizer, pressure cooker; operates between 116-126 <sup>0</sup> C; 10-20 psi; aluminum sterilizer has a 30 x 32.2 cm chamber; is supplied with chamber, lid with pressure gauge, immersion heater and safety valve, electric	NA	500.00	Sterilizing media and instrument
1 ea	Sterilizer, Autoclave; operates at 121 <sup>0</sup> C with dial for fast or slow exhaust; 0-60 minute timer; stainless unit with 66 41 cm chamber	NA	6500.00	Sterilizing media and instruments
1 ea	Sterilizer, dry heat with glass beads 120 V (S636) or 240 V (S637)	S636	395.00	Sterilizes instruments in hood between transfers
1 pkg	Stirring bars, magnetic; Teflon covered; leakproof, seamless and contain permanent Alnico V magnets	NA	58.00	Used for mixing stock solutions and in media preparation
1 ea	Stir bar retriever; contains an Alnico V magnet sealed in polyethylene; 18" length	NA	9.00	Retrieving stir bars from mixing vessel

## Items Required to Start a Plant Cell Culture Laboratory

QTY	ITEM DESCRIPTION	PHYTOTECH CATALOGUE NO.	APPROX. COST	ITEM FUNCTION
1 roll	Tape, autoclave indicator; impregnated to show the word "Autoclaved" after 15 minutes of exposure at 121 <sup>0</sup> C	T 998	3.95	Identifying autoclaved media
1 roll	Tape, Label; all purpose, self adhesive tape can be written on with pen or permanent marker			Labeling cultures, storage bottles, media vessels, etc.
2 ea	Thermometer; -20-150 <sup>0</sup> C temperature range	NA	40.00	Measuring temperature of liquids and culture room
1	Timer; electronic, countdown timer alarm, stopwatch feature	NA	30.00	Timing sterilization and general lab use
1 case	Towels; commercial, single fold		13.00	Can be sterilized to provide sterile work surface for cutting explants, gen lab use
500 ea	Weighing boats, plastic; can be utilized for liquid or solid samples; 1-5/8" x 1-5/8" x 5/16"	W 879	26.50	Measuring chemicals
500 ea	Weighing boats, plastic; can be utilized for liquid or solid samples; 3-5/16" x 3-5/16" x 1"	W 880	32.75	Measuring chemicals
500 ea	Weighing boats, plastic; can be utilized for liquid or solid samples; 5-1/2" x 5-1/2" x 1"	W 881	75.75	Measuring chemicals

The glassware used in tissue culture can generally be found in most laboratories. The glassware, particularly the culture vessels, should be made of Pyrex or borosilicate glass. Due to the increasing expense of this type of glass, many laboratories are successfully converting to soda glass, which may be seven to eight times cheaper. Wide-neck Erlenmeyer flasks (50-, 125-, 250-ml capacity) are commonly used as culture vessels; large volume Erlenmeyer flasks are required for media preparation. Test tubes, petri dishes, mason jars, baby food jars, and other glassware can also be adapted to tissue culture. Since all new glass may release substances that affect the composition of the medium, it is recommended that all new glassware be filled with water, autoclaved twice with detergent, washed, and rinsed between washes before being used for tissue culture. Other glassware commonly required in a tissue culture facility includes beakers, volumetric flasks, pipettes, and graduated cylinders.

